

# Package ‘SpectroX’

January 11, 2018

**Type** Package

**Title** Proteotypic peptide selection and SRM/PRM/HRM panel creation

**Version** 2.0

**Author** Erik Ahrne

**Maintainer** Erik Ahrne <erik.ahrne@unibas.ch>

**Description** Select proteotypic peptides based on full proteom analysis MaxQuant 1.6 search results or create SRM/PRM/HRM panel/library from MaxQuant 1.6 search results.

**License** GPL-3

**Encoding** UTF-8

**LazyData** true

**Suggests** testthat

**RoxygenNote** 6.0.1

## R topics documented:

barplotPeptidesPerProtein . . . . .	2
barplotPetideCountPerProtein . . . . .	2
createComplementaryIsotopeLibrary . . . . .	3
createLibrarySpectrum . . . . .	4
createSpectralLibrary . . . . .	4
digestProteome . . . . .	5
getCMDLineOptions . . . . .	6
getEmpiricalIRT . . . . .	7
getFragmentSequence . . . . .	8
getIRTModel . . . . .	8
getPeptides . . . . .	9
getPTMColors . . . . .	10
getSearchedModifications . . . . .	11
getUserOptions . . . . .	12
parseMaxQuantMSMS . . . . .	12
parseRange . . . . .	14
parseTargetsFile . . . . .	14
plotIRTCalibration . . . . .	15
proteotypicPeptideExport . . . . .	16
skylineExport . . . . .	17
spectroDiveExport . . . . .	17
spectronautExport . . . . .	18

**Index****19**

---

`barplotPeptidesPerProtein`*plot ranking metric vs peptide (per protein) barplot*

---

**Description**

plot ranking metric vs peptide (per protein) barplot

**Usage**

```
barplotPeptidesPerProtein(df, pepLenTrunc = 12, pepLabCex = 0.7,  
  rankingMetric = "rankingMetric", ptmCol = getPTMColors(levels(df$ptm)),  
  ...)
```

**Arguments**

<code>df</code>	data.frame
<code>pepLenTrunc</code>	integer AFADAMEVIPSTLAENAGLNPISTVTELRL -> AFADAMEVIPSTLAE..
<code>pepLabCex</code>	default 0.7
<code>rankingMetric</code>	character
<code>ptmCol</code>	data.frame row.names: ptm col:col

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

`barplotPetideCountPerProtein`*plot peptide count per protein*

---

**Description**

plot peptide count per protein

**Usage**

```
barplotPetideCountPerProtein(spectralLibrary, protLabCex = 0.9,  
  acLenTrunc = 12, col = "blue", ...)
```

**Arguments**

spectralLibrary	data.frame
protLabCex	default 0.9
acLenTrunc	integer default 12 "SOMEVERYLONGAC" -> "SOMEVERY.."

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

```
createComplementaryIsotopeLibrary
```

*Create complementary isotope (Arg, Lys H/L) spectral library*

---

**Description**

Create complementary isotope (Arg, Lys H/L) spectral library

**Usage**

```
createComplementaryIsotopeLibrary(sl)
```

**Arguments**

sl	tibble or data frame
----	----------------------

**Details**

No details

**Value**

data.frame

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

createLibrarySpectrum	<i>Parse spectrum framgment match information listed in MaxQuant 'Identifications' results</i>
-----------------------	------------------------------------------------------------------------------------------------

---

**Description**

Parse spectrum framgment match information listed in MaxQuant 'Identifications' results

**Usage**

```
createLibrarySpectrum(psm)
```

**Arguments**

psm	a row of MaxQuant 'Identifications' data.frame
-----	------------------------------------------------

**Details**

No details

**Value**

list with attributes intensity, ionType, charge, mz, isNL

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

createSpectralLibrary	<i>Create Spectral Library</i>
-----------------------	--------------------------------

---

**Description**

Create Spectral Library

**Usage**

```
createSpectralLibrary(tb, minFragNb = 3, minNbTransitions = 5,  
  maxNbTransitions = 5, minBasePeakFraction = 0, ionTypeFilter = c("a",  
  "b", "x", "y"), includeNL = T, rankingMetric = "adjustedIntensitySum")
```

**Arguments**

tb	tibble maxQuant spectrum level search results
minFragNb	min frag number default 3 (i.e. b3 and y3 will be kept)
minNbTransitions	minimum number of fragments
minBasePeakFraction	minimum intensity fraction of base peak (most intense kept fragment)
ionTypeFilter	selected ion type default a,b,x,y
includeNL	include neutral loss peaks default TRUE
rankingMetric	character. column name, used for ranking peptide, default "adjustedIntensitySum", c("adjustedIntensitySum", "psmScore", "precApexIntensity")

**Details**

No details

**Value**

data.frame

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

digestProteome	<i>Get peptide candidates per protein</i>
----------------	-------------------------------------------

---

**Description**

Get peptide candidates per protein

**Usage**

```
digestProteome(proteins, proteaseRegExp = "[KR](?!P)", dispProgressBar = T,  
  peptideLengthRange = c(6, 21), nbMiscleavages = 0,  
  exclusivePeptides = T, trimAC = T)
```

**Arguments**

proteins list  
 proteaseRegExp character '[KR](?!P)' - trypsin  
 dispProgressBar T/F  
 peptideLengthRange integer vector default c(6,21)  
 nbMisleavages 0  
 exclusivePeptides TRUE/FALSE get exclusive peptides only i.e. peptides mapping to a single protein. default T  
 trimAC TRUE/FALSE default TRUE splP62258|1433E\_HUMAN -> P62258

**Value**

data.frame peptide,protein

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

getCMDLineOptions	<i>Define and User Specified Command Line Options</i>
-------------------	-------------------------------------------------------

---

**Description**

Define and User Specified Command Line Options

**Usage**

```
getCMDLineOptions(version = version)
```

**Arguments**

version SpectroX version number

**Details**

No details

**Value**

list() cmd line options

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

getEmpiricalIRT	<i>Add iRT metric to data.frame</i>
-----------------	-------------------------------------

---

**Description**

Add iRT metric to data.frame

**Usage**

```
getEmpiricalIRT(tb, fit)
```

**Arguments**

tb	tibble containing colum labelled "Retention time"
fit	lm object

**Details**

No details

**Value**

vector of normalised rt (empirical irt)

**Note**

No note

**References**

Using iRT, a normalized retention time for more targeted measurement of peptides, Escher et al. (2012), <http://www.ncbi.nlm.nih.gov/pubmed/22577012>

**Examples**

```
print("No examples")
```

---

getFragmentSequence	<i>Get Fragment Sub Sequence</i>
---------------------	----------------------------------

---

**Description**

Get Fragment Sub Sequence

**Usage**

```
getFragmentSequence(peptide = peptide, ionType = ionType,  
  fragmentNb = fragmentNb)
```

**Arguments**

peptide	character string
ionType	c("a","b","x","y")
fragmentNb	integer

**Details**

No details

**Value**

character string

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

getIRTModel	<i>Get linear model predicting iRT as a function retention time (column name "rt")</i>
-------------	----------------------------------------------------------------------------------------

---

**Description**

Get linear model predicting iRT as a function retention time (column name "rt")

**Usage**

```
getIRTModel(tb)
```



**Arguments**

tb                      tibble including column "Retention time" and "Sequence"

**Details**

Uses Robust Linear Regression. Retention times of at least 3 iRT unique peptides needs to be listed in the input tibble

**Value**

list including fit rlm object, data.frame peptide,rt,irtRef

**Note**

No note

**References**

Using iRT, a normalized retention time for more targeted measurement of peptides, Escher et al. (2012), <http://www.ncbi.nlm.nih.gov/pubmed/22577012>

**Examples**

```
print("No examples")
```

---

getPeptides

*Get Mz Shift of complementary peptide/fragment ion*


---

**Description**

Get Mz Shift of complementary peptide/fragment ion

**Usage**

```
getPeptides(proteinSeq, proteaseRegExp = "[KR](?!P)", nbMiscleavages = 0)
```

**Arguments**

proteinSeq	protein sequence
proteaseRegExp	protease Regular Expression
nbMiscleavages	default 0
aaSeq	character string
charge	default 1
isHeavy	c(T,F)
annotSpectrum	Annotated Spectrum object

**Details**

Calculate mass shift of complimentary spectrum.

- PETIDEK (light) -> 8.014199
- PETIDEK (heavy) -> -8.014199
- PETIDER (light) -> 10.008269
- PETIDER (heavy) -> -10.008269

No details

No details

**Value**

numeric

annotSpectrum Annotated Spectrum object

vector of peptides

**Note**

No note

No note

No note

**References**

NA

NA

**Examples**

```
print("No examples")
Created complimentary Heavy/Light annotated spectrum with updates precursor and fragment m/z values.
print("No examples")
Digest protein
print("No examples")
```

---

getPTMColors

*assign colors to ptms*

---

**Description**

assign colors to ptms

**Usage**

```
getPTMColors(ptms)
```

**Arguments**

ptms

vector

**Value**

data.frame

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

getSearchedModifications

*Get list of variable modifications considered by MaxQuant*

---

**Description**

Get list of variable modifications considered by MaxQuant

**Usage**

```
getSearchedModifications(tb, ignoreArgLysIsoLabel = T)
```

**Arguments**

ignoreArgLysIsoLabel	default T Do not consider Arg Lys C,N heavy isotope label
df	tibble or data.frame

**Details**

No details

**Value**

character vector of modification names (does not return "Arg10", "Lys8" labels )

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

getUserOptions	<i>check and return User Specified Command Line Options</i>
----------------	-------------------------------------------------------------

---

**Description**

check and return User Specified Command Line Options

**Usage**

```
getUserOptions(cmdlineOptions)
```

**Arguments**

list()	command line options
--------	----------------------

**Details**

No details

**Value**

list() user options

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

parseMaxQuantMSMS	<i>Parse MaxQuant msms.txt</i>
-------------------	--------------------------------

---

**Description**

Parse MaxQuant msms.txt

**Usage**

```
parseMaxQuantMSMS(file, pepCutoff = 0.05, targetPeptides = NA,  
  targetProteins = NA, filterContaminants = T,  
  contaminantRegExp = "^CON_", selectedPTMRegExp = NA,  
  filterNonExclusivePeptides = T, pepLength = c(1, Inf),  
  chargeState = c(1, Inf), label = NA, maxMissedCleavages = 0,  
  keepBestSpectrumOnly = T, requiredPepSeqRegExp = ".", ...)
```

## Arguments

<code>file</code>	<code>path</code>
<code>pepCutoff</code>	numeric default 0.05
<code>targetPeptides</code>	character default NA
<code>targetProteins</code>	character default NA
<code>filterContaminants</code>	TRUE
<code>contaminantRegExp</code>	'^CON_'
<code>selectedPTMRegExp</code>	NA
<code>filterNonExclusivePeptides</code>	default TRUE
<code>pepLength</code>	default [1,Inf)
<code>chargeState</code>	default [1,Inf)
<code>label</code>	(Arg10 and Lys8 only filter) options NA (deafault,no filter), 'light','heavy'
<code>keepBestSpectrumOnly</code>	keep top-scoring spectrum only, per peptidoform, default TRUE
<code>requiredSequenceRegExp</code>	dafault '.' (no filter)

## Details

No details

## Value

data.frame of maxQuant psm level results

## Note

No note

## References

NA

## Examples

```
print("No examples")
```

---

parseRange	<i>Parse range in format 2:4 and return integer vector. Inthis case 2,3,4</i>
------------	-------------------------------------------------------------------------------

---

**Description**

Parse range in format 2:4 and return integer vector. Inthis case 2,3,4

**Usage**

```
parseRange(rangeStr)
```

**Arguments**

rangeStr            e.g. '2:4'

**Details**

No details

**Value**

vectot of integers

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

parseTargetsFile	<i>Parse list of target peptides/proteins</i>
------------------	-----------------------------------------------

---

**Description**

Parse list of target peptides/proteins

**Usage**

```
parseTargetsFile(file)
```

**Arguments**

file                    path

**Details**

No details

**Value**

list() peptides, proteins

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

plotIRTCalibration	<i>Plot IRT calibration curve</i>
--------------------	-----------------------------------

---

**Description**

Plot IRT calibration curve

**Usage**

```
plotIRTCalibration(irtModel, cex = 1.5, col = "blue", pch = 19, ...)
```

**Arguments**

irtModel	list
----------	------

**Details**

No details

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

```
proteotypicPeptideExport
```

```
# Export top X peptide per protein (ordered by 'adjustedIntensitySum')  
and add theoretical peptides ranked by length if fewer than X peptides  
were identified by MaxQuant
```

---

## Description

# Export top X peptide per protein (ordered by 'adjustedIntensitySum') and add theoretical peptides ranked by length if fewer than X peptides were identified by MaxQuant

## Usage

```
proteotypicPeptideExport(spectralLibrary, theoPeptides = NA,  
  targetProteins = unique(spectralLibrary$protein),  
  nbPeptidesPerProtein = 5, outFile = paste0(tempdir(), "/tmp.xls"))
```

## Arguments

spectralLibrary	tibble
theoPeptides	data.frame 'peptide' 'protein' 'length'
targetProteins	list of protein accession numbers e.g. P62258. default all proteins in spectralLibrary
nbPeptidesPerProtein	number of peptides per protein to be exported
outFile	path to output file default tmp.xls in tmp dir

## Note

No note

## References

NA

## Examples

```
print("No examples")
```



---

skylineExport	<i>Write skyline compatible xls file</i>
---------------	------------------------------------------

---

**Description**

Write skyline compatible xls file

**Usage**

```
skylineExport(sl, file)
```

**Arguments**

sl                      tibble or data frame

**Details**

No details

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

spectroDiveExport	<i>Write SpectroDive compatible xls file</i>
-------------------	----------------------------------------------

---

**Description**

Write SpectroDive compatible xls file

**Usage**

```
spectroDiveExport(sl, file)
```

**Arguments**

sl                      tibble or data frame

**Details**

No details

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

spectronautExport	<i>Write Spectronaut compatible xls file</i>
-------------------	----------------------------------------------

---

**Description**

Write Spectronaut compatible xls file

**Usage**

```
spectronautExport(sl, file)
```

**Arguments**

sl	tibble or data frame
----	----------------------

**Details**

No details

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

# Index

barplotPeptidesPerProtein, [2](#)  
barplotPetideCountPerProtein, [2](#)  
  
createComplementaryIsotopeLibrary, [3](#)  
createLibrarySpectrum, [4](#)  
createSpectralLibrary, [4](#)  
  
digestProteome, [5](#)  
  
getCMDLineOptions, [6](#)  
getEmpiricalIRT, [7](#)  
getFragmentSequence, [8](#)  
getIRTModel, [8](#)  
getPeptides, [9](#)  
getPTMColors, [10](#)  
getSearchedModifications, [11](#)  
getUserOptions, [12](#)  
  
parseMaxQuantMSMS, [12](#)  
parseRange, [14](#)  
parseTargetsFile, [14](#)  
plotIRTCalibration, [15](#)  
proteotypicPeptideExport, [16](#)  
  
skylineExport, [17](#)  
spectroDiveExport, [17](#)  
spectronautExport, [18](#)