

# Package ‘SpectroX’

December 6, 2017

**Type** Package

**Title** Proteotypic peptide selection and SRM/PRM/HRM panel creation

**Version** 1.0

**Author** Erik Ahrne

**Maintainer** Erik Ahrne <erik.ahrne@unibas.ch>

**Description** Select proteotypic peptides based on full proteom analysis MaxQuant 1.6 search results or create SRM/PRM/HRM panel/library from MaxQuant 1.6 search results.

**License** GPL-3

**Encoding** UTF-8

**LazyData** true

**Suggests** testthat

**RoxygenNote** 6.0.1

## R topics documented:

createAnnotatedSpectrum . . . . .	1
getEmpiricalIRT . . . . .	2
getFragmentSequence . . . . .	3
getIRTModel . . . . .	4
getLabelMzShift . . . . .	4
getSearchedModifications . . . . .	5
parseMaxQuantMSMS . . . . .	6
plotIRTCalibration . . . . .	7

<b>Index</b>	<b>8</b>
--------------	----------

---

createAnnotatedSpectrum

*Parse spectrum framgment match information listed in MaxQuant  
'Identifications' results*

---

## Description

Parse spectrum framgment match information listed in MaxQuant 'Identifications' results

**Usage**

```
createAnnotatedSpectrum(psm)
```

**Arguments**

psm                      a row of MaxQuant 'Identifications' data.frame

**Details**

No details

**Value**

list with attributes intensity, ionType, charge, mz, isNL

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

getEmpiricalIRT	<i>Add iRT metric to data.frame</i>
-----------------	-------------------------------------

---

**Description**

Add iRT metric to data.frame

**Usage**

```
getEmpiricalIRT(tb, fit)
```

**Arguments**

tb                      tibble containing colum labelled "Retention time"  
fit                      lm object

**Details**

No details

**Value**

vector of normalised rt (empirical irt)

**Note**

No note

**References**

Using iRT, a normalized retention time for more targeted measurement of peptides, Escher et al. (2012), <http://www.ncbi.nlm.nih.gov/pubmed/22577012>

**Examples**

```
print("No examples")
```

---

getFragmentSequence	<i>Get Fragment Sub Sequence</i>
---------------------	----------------------------------

---

**Description**

Get Fragment Sub Sequence

**Usage**

```
getFragmentSequence(peptide = peptide, ionType = ionType,  
                    fragmentNb = fragmentNb)
```

**Arguments**

peptide	character string
ionType	c("a","b","x","y")
fragmentNb	integer

**Details**

No details

**Value**

character string

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

getIRTModel	<i>Get linear model predicting iRT as a function retention time (column name "rt")</i>
-------------	--

---

**Description**

Get linear model predicting iRT as a function retention time (column name "rt")

**Usage**

```
getIRTModel(tb)
```

**Arguments**

tb                      tibble including column "Retention time" and "Sequence"

**Details**

Uses Robust Linear Regression. Retention times of at least 3 iRT unique peptides needs to be listed in the input tibble

**Value**

list including fit rlm object, data.frame peptide,rt,irtRef

**Note**

No note

**References**

Using iRT, a normalized retention time for more targeted measurement of peptides, Escher et al. (2012), <http://www.ncbi.nlm.nih.gov/pubmed/22577012>

**Examples**

```
print("No examples")
```

---

getLabelMzShift	<i>Get Mz Shift of complementary peptide/fragment ion</i>
-----------------	---

---

**Description**

Get Mz Shift of complementary peptide/fragment ion

**Usage**

```
getLabelMzShift(aaSeq, charge = 1, isHeavy = T)
```

**Arguments**

aaSeq	character string
charge	default 1
isHeavy	c(T,F)

**Details**

Calculate mass shift of complimentary spectrum.

- PETIDEK (light) -> 8.014199
- PETIDEK (heavy) -> -8.014199
- PETIDER (light) -> 10.008269
- PETIDER (heavy) -> -10.008269

**Value**

numeric

**Note**

No note

**References**

NA

**Examples**

```
print("No examples")
```

---

`getSearchedModifications`*Get list of variable modifications considered by MaxQuant*

---

**Description**

Get list of variable modifications considered by MaxQuant

**Usage**

```
getSearchedModifications(tb)
```

**Arguments**

df	tibble or data.frame
----	----------------------

**Details**

No details

Value

character vector of modification names

Note

No note

References

NA

Examples

```
print("No examples")
```

---

parseMaxQuantMSMS	<i>Parse MaxQuant msms.txt</i>
-------------------	--------------------------------

---

Description

Parse MaxQuant msms.txt

Usage

```
parseMaxQuantMSMS(file, pepThrs = 0.05, targetPeptides = NA,
  targetProteins = NA, filterContaminants = T,
  contaminantRegExp = "^CON_", selectedPTMRegExp = NA,
  filterNonExclusivePeptides = T, minPepLength = 0, chargeState = 1:10,
  label = NA, maxMissedCleavages = 0)
```

Arguments

file	path
pepThrs	numeric default 0.05
targetPeptides	character default NA
targetProteins	character default NA
filterContaminants	TRUE
contaminantRegExp	'^CON_'
selectedPTMRegExp	NA
filterNonExclusivePeptides	default TRUE
minPepLength	default 0
chargeState	default [1,10]
label	(Arg10 and Lys8 only filter) options NA (deafault,no filter), 'light','heavy'

## Details

No details

## Value

data.frame of maxQuant psm level results

## Note

No note

## References

NA

## Examples

```
print("No examples")
```

---

plotIRTCalibration	<i>Plot IRT calibration curve</i>
--------------------	-----------------------------------

---

## Description

Plot IRT calibration curve

## Usage

```
plotIRTCalibration(irtModel, cex = 1.5, col = "blue", pch = 19, ...)
```

## Arguments

irtModel	list
----------	------

## Details

No details

## Note

No note

## References

NA

## Examples

```
print("No examples")
```

# Index

`createAnnotatedSpectrum`, [1](#)

`getEmpiricalIRT`, [2](#)

`getFragmentSequence`, [3](#)

`getIRTModel`, [4](#)

`getLabelMzShift`, [4](#)

`getSearchedModifications`, [5](#)

`parseMaxQuantMSMS`, [6](#)

`plotIRTCalibration`, [7](#)